


Multiplex qPCR and Cannabis Microbiome Sequencing Reveals Several Bacteria and Fungi Native to Cannabis Flowers

Kevin McKernan, Jessica Spangler, Lei Zhang, Vasisht Tadigotla, Yvonne Helbert, Kellie Dodd, Douglas Smith Medicinal Genomics Corporation

The Center for Disease Control estimates 128,000 people in the U.S. are hospitalized annually due to food-borne illnesses. As a result, the detection of mold and bacteria on agricultural products has become an important safety consideration. This risk extends itself to medical Cannabis and is of particular concern with inhaled, vaporized and even concentrated Cannabis products. As a result, third-party microbial testing has become a regulatory requirement in the medical and recreational Cannabis markets. Medicinal Genomics has developed a novel PCR-based assay for the detection of pathogenic microbes in Cannahis materials. This process consists of a proprietary system for DNA extraction called SenSATIVAx™ and a novel PCR assay called PathogINDICAtor™

PathogINDICAtor™ utilizes a novel qPCR-based assay that is contamination-free and provides an internal plant DNA control for every reaction. DNA detection is based on a 5' nuclease assay that directly measures the amount of plant and microbe DNA in a given sample. This technique provides robust sensitivity, specificity, and multiplexing capability. Positive results can then be confirmed through microbiome sequencing. PathogINDICAtor™ is able to provide a contamination-free technique through a proprietary process known as DREAM PCR1,2 which uses methylation-specific restriction enzymes to prevent the carryover of amplified products from one reaction to the next. We present



MGC Offerings resence/Absence Tests Threshold Tests Total Yeast & Mold (18S) E.coli STEC only Examples of PathogINDICAtor™ Data

Figure 4: Sample with no microbe DNA Figure 3: Sample containing plant and microbe DNA (Plant DNA = Green, Microbial DNA = esent. (Plant DNA = Green, Microbial blue). X-axis is cycles (every 1.5 minutes). Y-axis

Example 1: Comparison of a Visually-Contaminated Sample



s relative fluorescence units in a log 10 scale.





Figure 5: With only 1% of microbes able to grow into colonies³, there is the potential for false negative results with conventional culture-based methods. The cannabis sample has visible mold growing, however, the culture-based method does not detect the microbes.

PathogINDICAtorTM detects and fails the sample.

Example 2: Comparison of Multiple Methods








Figure 6: To aid in validation of SenSATIVAx™ and PathoglNDICAtor™, we utilized three alternative culture-based assays (3M™Petrifilm™, SimPlate®, and Biolumix). When a sample does not correlate we reflex to next-generation sequencing of

The Microbiome Profile of the Sequenced Positive Result

Varme		(% of classified reads)	411
	Byssochiamys spectabilis No. 5	18859 (48.49%)	Readcount Taxonomic name 0 Distarrota
	Eurotiales	18021 (46.34%)	0 - Opistyskonta
	Wolfigoria cocos MD-104 SS10	509 (1.31%)	12 — Fungi
	saccharomyceta	360 (0.93%)	358 — Dionya 6 — Ascernosta
	Disarva	358 (0.92%)	360 — sacchiromyceta
	leotiomyceta	286 (0.74%)	1 Pezieomycotina
			296 lectiomyceta
	Pseudogymnoascus destructans 20631-21	97 (0.25%)	0 Eurotiomycetes
	Penicillium	72 (0.19%)	0 Eurodomycetidae 19021 Furodoles
	Saccharomycetales	43 (0.11%)	0 Thermoscaceae
	Neotyphodium aotearoae	22 (0.05%)	0 Byssochlarnys
			0 Byssochlamys spi
	(Remaining organisms)	263 (0.68%)	18859 Byssochlamys s

fiseq with our Yeast & Mold Detection Assay to reveal the microbes contained in the sample. Using OneCodex to analyze the data, this sample is primarily Byssochlamys, which produces mycotoxin and should be detected for safety. (Sample 3 in Figure 8 is also a representation of this sample) In general, this data can be used for many purpose, most notably to determine the possible contaminates in a greenhouse



Figure 8: A look at multiple sequenced PathogINDICAtor™ positive cannabis samples and the microbiomes found on samples and the microbiomes found on each. This data can help in determining trends within multiple cannabis samples. Notice that 9 of 10 have Penicillium as a part of the microbiome, which is a known

encopnyte".

Sample 1 is the Australian Bastard sample referenced in the presentation titled "Genomic, Terpene, and Cannabinoid Profiles of a Putatively Novel Cannabis

As many regulations are beginning to mandate "heat killing" of microbial content, we must remain aware of how these drying techniques often confound culturebased methods used to monitor colony forming units (CFU). Even though this "heat kill" process may be effective at sterilizing some of the culturable microbial content it does not eliminate various pathogenic toxins like Aflatoxin or the DNA that encodes the Aflatoxin gene. This is of particular importance as Aflatoxin is a carcinogen. The clearance of Aflatoxin requires the liver enzyme CYP3A4 and this liver enzyme is potently inhibited by Cannabinoids⁵⁻⁷. With the publication of the Cannabis genome and many pathogenic microbial genomes, we designed several multiplexed quantitative PCR (qPCR) assays to detect pathogenic DNA in a background of host Cannabis DNA.

One of these qPCR assays is designed to detect total yeast and mold and thus targets the 18S rDNA ITS (Internal Transcribed Spacer) regions. ITS regions are routinely used to sequence the microbiomes of samples and identify and itemize the collection of microbial communities present in a given sample. The addition of next-generation sequencing primer tails to this yeast and mold qPCR assay enables reflexive sequencing of samples testing positive for yeast and mold. Nextgeneration sequencing thus reveals a digital record of the microbiological community on a given Cannabis sample. Cannabis samples that tested positive for yeast and mold with gPCR and culture-based techniques were sequenced to precisely identify the collection of microbes present in the failing gPCR samples Microbes harmful to both the plant and humans were identified (Botrytis, Magneporthe, Fusarium and Aspergillus) while over 180 different microbes including endophytes were also present. This qPCR assay presents a real-time tool for both quality testing and geospatial monitoring of microbial communities in cannabis production facilities. This highly sensitive assay allows for the detection of very low levels of contamination. This can be used to proactively monitor for outbreaks before they can spread.

Research and Development Contact Business Development Contact

Email: Jessica.Spangler@medicinalgenomics.com

Kellie Dodd

Email: Kellie.Dodd@medicinalgenomics.com

Sales Contact

John Geanaconoulos Email: ig@medicinalgenomics.com

References



www.medicinalgenomics.com