

# Paw-Pads™

NON-SLIP TRACTION PADS

# Paw-Pads™

NON-SLIP TRACTION PADS

# Paw-Pads™

NON-SLIP TRACTION PADS

**Paw-Pads™ Non-Slip Traction Pads** helps senior and arthritic, dogs from injury due to slipping on hardwood, tile, and other slippery surfaces.

Dogs prefer our traction pads over dog boots dog socks, dog rubber boots and dog shoes, because unlike dog boots, dog shoes and dog rubber socks, your dog won't even know they have paw pads on as it provides traction for senior and arthritic dogs and helps them navigate surfaces that are slippery.

**Paw Pads™** are made from medical grade material. The ultra-thin fabric allows paws to breathe and provides excellent traction on tile, hardwood, marble, vinyl, or any slippery surfaces.



No more slipping on hardwood floors. Paw-Pads™ also provides traction for pre op and post op dogs, dogs with hip dysplasia, senior dogs with arthritis, and dogs recuperating from hip or back surgery.



Questions

[sales@pupgearpetproducts.com](mailto:sales@pupgearpetproducts.com)

## New and Improved Sleeker-Slimmer-Better Traction



The Next Generation of Paw Pads™ is 30% thinner with 25% better traction and 25% stronger adhesive. Our New Material and Improved Adhesive make a difference in the quality of life for senior dogs and dogs recuperating for ACL or Hip surgery.

Paw Pads™ are made with 3M® Skin Safe Adhesive that is Non-Toxic and Hypo- Allergenic.

Paw Pads™ are so comfortable dogs do not realize they have them on.

Perforations allow the pad to breath while providing traction on slick surfaces like hardwood and tile.

Paw Pads™ are similar to Band-Aids and the length of time they remain on a dogs pad depends on their environment. Typically, they will stay on for 2-7 days.

Paw Pads™ can be reapplied using our adhesive kit or medical grade super glue.

Veterinarian Approved



LIFESTYLE PRODUCTS APPROVED BY DOGS

1.877.781.PUPS

# Paw Pads™

NON-SLIP TRACTION PADS

# Paw Pads™


NON-SLIP TRACTION PADS

# Paw Pads™

NON-SLIP TRACTION PADS

## Paw-Pads™ Measuring Guide.

If your dog's paws are between sizes indicated below, order the smaller size of Paw-Pads™.



Paw-Pads™ do not need to cover the entire pad to be effective. It may also be necessary to trim the paw-pad to get proper fit.

### Choosing the right size for you dogs paws.

Measure your dogs paw from the top of the digital pad to the bottom of the metatarsal pad. Pick the size that is smaller than your dogs overall paw size.

If the Length of a dogs paws is between 3/4"-7/8" Use SM/Medium

If the Length of a dogs paws is between 7/8"-2" Use Large

## Important Tips

It is not necessary to apply pads to all the digits on all your dogs paws. Usually two or three pads per foot will provide good traction.

Paw Pads do not need to cover the entire pad to be effective, it maybe necessary to trim the paw pad to get proper fit.

The condition of a dogs pad may affect the adhesion. Smooth paws will adhere better than rough or cracked pads.

Paw Pads can be reapplied using the Paw Pad Adhesive Kit that includes grade medical grade super glue These adhesives are non-toxic and are used for bonding skin together in place of stitches and staples and perfectly safe to use with Paw Pads™

**DO NOT USE Paw Pads™ on open wounds unless directed by veterinarian. Do not apply to pads that are irritated or swollen. Check your dog's paws for signs of adverse reaction, allergy or skin related conditions.**



Apply one pad at a time. Remove from adhesive strip and apply directly to dogs pad. Paw Pads are pressure sensitive and pressure must be applied to the Paw Pad to achieve proper adhesion. Pressure can be applied by holding the dogs paw to your hand or having your dog stand without moving for at least 90 seconds.

## Paw Pad™ Adhesive Kit

Get More Mileage out of Paw Pads™. Use to re-secure Paw Pads™ if needed.

Paw Pad™ Adhesive Kit is a cost effective means to extend the life of Paw Pads. Kit includes everything you need to make reapplying easy! The Kit includes:

- 3 ml Medical Grade Veterinary Skin Adhesive
- 5' Plastic Forceps
- 2- 1/2" Plastic Applicators



Questions  
sales@pupgearpetproducts.com

1.877.781.PUPS



LIFESTYLE PRODUCTS APPROVED BY DOGS